

# ORTOVIA

Upper Extremity

ORTOVIA

*"Smart Fixation"*



Dorsal-side  
fixing



radiolucent  
Joint



telescopic  
Body



Sliding  
clamp



Easy  
lengthening



MRI  
Compatibility

## Fixat Types

- Wrist
- Forearm
- Humerus


ORTOVIA Upper Extremity Fixator System is a series of modular monolateral external fixators for upper extremity that was developed under many years of clinical experience. The system is composed of wrist, forearm and humerus fixators.

## Features

- Modularity
- Easy application
- Short surgery time
- MRI compatibility
- Flexibility-versatility
- Stability & safety in corrections
- Patient's comfort
- Radiolucent components
- Short learning curve


## Indications

- Bone loss
- Open fractures
- Comminuted fractures
- Non-union
- Angular deformities
- Limb reconstruction




# Specifications

ORTOVIA

Reference No.	Description	Image
8680628122357	Ortovia Sliding Clamp Small	
8680628122371	Ortovia Metaphysis Clamp Small	
8680628122388	Ortovia Polyclamp Small	
8680628122401	Ortovia Open Clamp Small	
8680628136439	Ortovia 3D Joint Small	
8680628136477	Ortovia Hinge Joint Small	
8680628136507	Ortovia Fixator Body LS10	
8680628136514	Ortovia Fixator Body LS20	
8680628136521	Ortovia Fixator Body LS25	
8680628122999	Ortovia Wrist Fixator Short	
8680628123002	Ortovia Wrist Fixator Long	

Reference No.	Description	Image
8680628123705	Ortovia Shanz Screw 20mm Thread (S.Steel) Ø3X50MM	
8680628123712	Ortovia Shanz Screw 20mm Thread (S.Steel) Ø3X60MM	
8680628123729	Ortovia Shanz Screw 20mm Thread (S.Steel) Ø3X70MM	
8680628123033	Ortovia Shanz Screw 35mm Thread (S.Steel) Ø3X70MM	
8680628123040	Ortovia Shanz Screw 35mm Thread (S.Steel) Ø4X70MM	
8680628123057	Ortovia Shanz Screw 35mm Thread (S.Steel) Ø4X80MM	
8680628124405	Ortovia Shanz Screw 20mm Thread (Titanium) Ø3X50MM	
8680628124412	Ortovia Shanz Screw 20mm Thread (Titanium) Ø3X60MM	
8680628124429	Ortovia Shanz Screw 20mm Thread (Titanium) Ø3X70MM	
8680628123736	Ortovia Shanz Screw 35mm Thread (Titanium) Ø3X70MM	
8680628123743	Ortovia Shanz Screw 35mm Thread (Titanium) Ø4X70MM	
8680628123750	Ortovia Shanz Screw 35mm Thread (Titanium) Ø4X80MM	



Headquarter & Sales:  
Metroport Buidence No :143, E5 Üzeri Kültür  
Sokak No: 1, 34180 Bahçelievler / Istanbul  
Tel: +90 212 442 50 00 Fax: +90 212 442 46 34  
info@spinamer.com, sales@spinamer.com  
www.spinamer.com

R&D Department:  
Yıldız Technical University  
Davutpasa Campus, Technology Development Zone  
B2 Block, No: 114, Esenler / Istanbul  
Tel: +90 212 442 50 00 Fax: +90 212 442 46 34  
E-mail: research@spinamer.com

Factory:  
Sanayi Mah. 1715, Sok. No: 32  
34538 Esenyurt / Istanbul  
Tel: +90 212 44 16 444 Fax: +90 212 442 46 34  
E-mail: factory@spinamer.com

