

Ortovia android

ORTOVIA

External Lengthening System



Software Control



Remote Access



Internet Connection



Lengthening



Shortening



Precise Mechatronic

Ortovia android Lengthening system designed to be used in lengthening and shortening purpose of the limb is a rigid and precision mechatronic external lengthening system. Ortovia app. is a special android-based software that authorized orthopedic surgeon can program and operate the system with remote access via the Internet.


On this system Ortovia.app is installed in patient's mobile phone, the doctor by his mobile phone or PC can reach the relevant patient's lengthening system and obtain information about the current situation, and also he can reprogram the lengthening or shortening process.

Indications:

- lengthening and shortening process in long bones


Contraindications:

- In case the patient does not have enough psychological and intellectual abilities to sustain the operation and postoperative treatment .
- Other known contraindications of limb lengthening process.



Specifications

ORTOVIA

Reference No.	Description	Image
100.8030	Ortovia L30 android fixator	
100.8025	Ortovia L25 android fixator	
100.8020	Ortovia L20 android fixator	
100.8125	Ortovia LS25 android fixator	
100.8120	Ortovia LS20 android fixator	
100.8225	Ortovia LR25 android fixator	
100.8220	Ortovia LR20 android fixator	
100.8320	Ortovia LRS20 android fixator	
100.4110	ORTOVIA SLYDING CLAMP	
100.4111	ORTOVIA SLYDING CLAMP 1/2	
100.4112	ORTOVIA METAPHYSIS CLAMP	
100.4113	ORTOVIA POLYCLAMP	
100.8800	Electronic Control Unit	
100.8900	Ortovia app. Software	



Headquarter & Sales:
Metroport Buidence No :143, E5 Üzeri Kültür
Sokak No: 1, 34180 Bahçelievler / Istanbul
Tel: +90 212 442 50 00 Fax: +90 212 442 46 34
info@spinamer.com, sales@spinamer.com
www.spinamer.com

R&D Department:
Yıldız Technical University
Davutpasa Campus, Technology Development Zone
B2 Block, No: 114, Esenler / Istanbul
Tel: +90 212 442 50 00 Fax: +90 212 442 46 34
E-mail: research@spinamer.com

Factory:
Sanayi Mah. 1715. Sok. No: 32
34538 Esenyurt / Istanbul
Tel: +90 212 44 16 444 Fax: +90 212 442 46 34
E-mail: factory@spinamer.com

